



Product Datasheet

Eph receptor A2/EPHA2 Antibody (orb402207)


Description	Eph receptor A2/EPHA2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, WB
Immunogen	E. coli-derived human Eph receptor A2 recombinant protein (Position: M851-N970).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	125 kDa
Uniprot ID	P29317
Dilution Range	Western blot, 0.1-0.5µg/ml Immunocytochemistry/Immunofluorescence, 5µg/ml Flow Cytometry, 1-3µg/1x10 ⁶ cells Direct ELISA, 0.1-0.5µg/ml
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of A549 cells us...



Flow Cytometry analysis of U20S cells us...



WB analysis of Eph receptor A2 using ant...