



Product Datasheet

ASXL1 Antibody (orb402231)

Description	ASXL1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthetic peptide corresponding to a sequence of human ASXL1 (KKERTWAEAAARLVLENYSAPMTKQILQVIEAE).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	165 kDa
Uniprot ID	Q8IXJ9
Dilution Range	Western blot, 0.1-0.5µg/ml Flow Cytometry, 1-3µg/1x10 ⁶ cells
Expiration Date	12 months from date of receipt.



Flow
Cytometry
analysis of
HepG2 cells
U...



WB analysis of
ASXL1 using
anti-ASXL1
an...