



Product Datasheet

Uteroglobin/Scgb1a1 Antibody (orb402269)


Description	Uteroglobin/Scgb1a1 Antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	E. coli-derived rat Uteroglobin recombinant protein (Position: S20-V96).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	10 kDa
Uniprot ID	P17559
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml ELISA(Cap), 1-5µg/ml
Expiration Date	12 months from date of receipt.



WB analysis of Uteroglobin using anti-Ut...



Sandwich ELISA - Recombinant rat SCGB1A1...



IHC analysis of Uteroglobin using anti-U...