

Product Datasheet

CTLA4 Antibody (orb402310)


Description

CTLA4 Antibody


Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	E. coli-derived mouse CTLA4 recombinant protein (Position: E36-D161). Mouse CTLA4 shares 67.5% amino acid (aa) sequence identity with human CTLA4.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	24993 MW
Uniprot ID	P09793
Dilution Range	Western blot, 0.1-0.5µg/ml, Mouse Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Mouse, Rat, By Heat ELISA , 0.1-0.5µg/ml, Mouse
Expiration Date	12 months from date of receipt.



WB analysis of CTLA4 using anti-CTLA4 an...



IHC analysis of CTLA4 using anti-CTLA4 a...



IHC analysis of CTLA4 using anti-CTLA4 a...