

Product Datasheet

IL33 Antibody (orb402322)

Description

IL33 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

E. coli-derived human IL33 recombinant protein (Position: A95-T270).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results.

By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users.

. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

18 kDa

Uniprot ID
095760
Dilution Range

Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml ELISA (Cap), 1-5µg/ml

Expiration Date

12 months from date of receipt.

100KD


55KD

35KD


25KD

15KD

WB analysis
of IL33 using
anti-IL33
anti...



Sandwich
ELISA -
Recombinant
human IL33 ...



IHC analysis
of IL33 using
anti-IL33
ant...