

Product Datasheet

DDX58 Antibody (orb402335)

Description

DDX58 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, ICC, IF, WB

Applications

Immunogen

E. coli-derived human DDX58 recombinant protein (Position: H871-K925). Human DDX58 shares 83.3% amino acid (aa) sequence identity with mouse DDX58.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

106600 MW

Uniprot ID


[O95786](#)

Dilution Range


Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat
Immunocytochemistry/Immunofluorescence, 2µg/ml, Human
Immunofluorescence, 5µg/ml, Human Flow Cytometry, 1-3µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



Flow Cytometry analysis of U87 cells us...



Flow Cytometry analysis of A431 cells us...



WB analysis of DDX58 using anti-DDX58 an...