

Product Datasheet

Fibroblast activation protein, alpha/FAP Antibody (orb402355)

Description Fibroblast activation protein, alpha/FAP Antibody

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, WB

Immunogen E. coli-derived human FAP recombinant protein (Position: R550-A619).

Form/Appearance Lyophilized

Concentration Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note For research use only

Application notes Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype Rabbit IgG


Clonality Polyclonal

MW 95 kDa

Uniprot ID **Q12884**

Dilution Range Western blot, 0.1-0.5µg/ml Direct ELISA, 0.1-0.5µg/ml

Expiration Date 12 months from date of receipt.



WB analysis of FAP using anti-FAP antibo...