

Product Datasheet

TRAIL/TNFSF10 Antibody (orb402364)

Description
TRAIL/TNFSF10 Antibody
Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

E. coli-derived human TRAIL recombinant protein
(Position: N40-A272).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration
of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After
reconstitution, at 4°C for one month. It can also be
aliquotted and stored frozen at -20°C for six months.
Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive
results. By Heat: Boiling the paraffin sections in 10mM
citrate buffer, pH6.0, for 20mins is required for the
staining of formalin/paraffin sections. Other applications
have not been tested. Optimal dilutions should be
determined by end users. Add 0.2ml of distilled water
will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


35 kDa

Uniprot ID
P50591
Dilution Range


Western blot, 0.1-0.5µg/ml, Human
Immunohistochemistry (Paraffin-embedded Section),
0.5-1µg/ml, Human, Mouse, Rat, By Heat

Expiration Date


12 months from date of receipt.



WB analysis
of TRAIL
using anti-
TRAIL an...



IHC analysis
of TRAIL
using anti-
TRAIL a...



IHC analysis
of TRAIL
using anti-
TRAIL a...