

Product Datasheet

PARK7/DJ1 Antibody (orb402399)

Description

PARK7/DJ1 Antibody

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

E. coli-derived mouse PARK7 / DJ1 recombinant protein (Position: A2-D189).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


22 kDa, 25 kDa

Uniprot ID
[Q99LX0](#)
Dilution Range


Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml ELISA (Cap), 1-5µg/ml

Expiration Date


12 months from date of receipt.



WB analysis of PARK7 / DJ1 using anti-PARK7 / DJ1 antibody.



IHC analysis of PARK7 / DJ1 using anti-PARK7 / DJ1 antibody.



IHC analysis of PARK7 / DJ1 using anti-PARK7 / DJ1 antibody.