

## Product Datasheet

### Periostin/Postn Antibody (orb402458)

## Description

## Periostin/Postn Antibody

### Species/Host

Rabbit

### Reactivity

Mouse, Rat

### Conjugation

Unconjugated

### Tested

ELISA, IHC, WB

### Applications

### Immunogen

E. coli-derived rat Periostin recombinant protein (Position: K577-Q810).

### Form/Appearance

Lyophilized

### Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

### Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

### Note

For research use only

### Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

### Isotype

Rabbit IgG

### Clonality

Polyclonal

### MW

90056 MW

### Uniprot ID

**D3ZAF5**


### Dilution Range

Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Direct ELISA, 0.1-0.5µg/ml


### Expiration Date

12 months from date of receipt.


1 2 3 4 5 6



WB analysis of Periostin using anti-Peri...



IHC analysis of Periostin using anti-Per...



IHC analysis of Periostin using anti-Per...