



Product Datasheet

Superoxide Dismutase 3/SOD3 Antibody (orb402479)


Description	Superoxide Dismutase 3/SOD3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence of human SOD3 (WTGEDSAEPNSDSA EWIRDMYAKVTEIWQEVMQRRDDD).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33 kDa
Uniprot ID	P08294
Dilution Range	Western blot, 0.25-0.5 µg/ml Immunohistochemistry(Paraffin-embedded Section), 0.5-1 µg/ml Immunocytochemistry/Immunofluorescence, 5 µg/ml
Expiration Date	12 months from date of receipt.



WB analysis of SOD3 using anti-SOD3 anti...



IF analysis of SOD3 using anti-SOD3 anti...



IHC analysis of SOD3 using anti-SOD3 ant...