

Product Datasheet

Neurocan/NCAN Antibody (orb402571)

Description

Neurocan/NCAN Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence of human Neurocan
(HENTWIGLNDRIVERDFQWTDNTGLQFENWRENQPDNFFA).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

200 kDa

Uniprot ID


[O14594](#)

Dilution Range


Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml

Expiration Date


12 months from date of receipt.



WB analysis of Neurocan using anti-Neuro...



IHC analysis of Neurocan using anti-Neur...



IHC analysis of Neurocan using anti-Neur...