

Product Datasheet

HOMER3 Antibody (orb402580)

Description

HOMER3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

E. coli-derived human HOMER3 recombinant protein (Position: R282-A360). Human HOMER3 shares 88.6% and 89.9% amino acid (aa) sequence identity with mouse and rat HOMER3, respectively.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


39836 MW

Uniprot ID
Q9NSC5
Dilution Range


Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By Heat Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat,

Expiration Date


12 months from date of receipt.



WB analysis of HOMER3 using anti-HOMER3 ...



IHC analysis of HOMER3 using anti-HOMER3...



IHC analysis of HOMER3 using anti-HOMER3...