

Product Datasheet

Epigen/Epgn Antibody (orb402590)

Description

Epigen/Epgn Antibody

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

E. coli-derived mouse Epigen recombinant protein (Position: A19-K109). Mouse Epigen shares 82.8% amino acid (aa) sequence identity with human Epigen.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

16799 MW

Uniprot ID


[Q924X1](#)

Dilution Range

ELISA , 0.1-0.5µg/ml, Mouse, Western blot, 0.1-0.5µg/ml, Rat, Mouse,

Expiration Date

12 months from date of receipt.



WB analysis of Epigen using anti-Epigen ...