


Product Datasheet

Rat BMP-9 ELISA Kit (orb402725)

Description	Rat BMP-9 ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: S320-R430
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	bone tissue, cell culture supernatants, serum and plasma (heparin, EDTA).
MW	47978 MW
Uniprot ID	M0RB87
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 120 pg/ml.
Expiration Date	Please enquire.



Rat BMP-9
ELISA Kit
standard
curve