


Product Datasheet

Human IL-23R/Interleukin-23 receptor ELISA Kit (orb402726)

Description	Human IL-23R/Interleukin-23 receptor ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: I25-I354
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q5VWK5
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 156 pg/ml.
Expiration Date	Please enquire.



Human IL-23R ELISA Kit standard curve