


Product Datasheet

Human Dystroglycan ELISA Kit (orb402728)

Description	Human Dystroglycan ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: H30-G653
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva and urine.
Uniprot ID	Q14118
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 500-1500 pg/ml.
Expiration Date	Please enquire.



Human
Dystroglycan
ELISA Kit
standard
cu...