


Product Datasheet

Mouse PGLYRP1/Pgrps ELISA Kit (orb402731)

Description	Mouse PGLYRP1/Pgrps ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: F19-E182
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	O88593
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is 33-100 ng/ml.
Expiration Date	Please enquire.



Mouse
PGLYRP1
ELISA Kit
standard
curve