


Product Datasheet

Mouse PON1/Paraoxonase 1 ELISA Kit (orb402734)

Description	Mouse PON1/Paraoxonase 1 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A30-L355
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
Uniprot ID	P52430
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:800, concentration in serum and plasma is 700-2000 ng/ml.
Expiration Date	Please enquire.



Mouse
PON1 ELISA
Kit standard
curve