



Product Datasheet

CD274 antibody (orb402781)

Description	Mouse monoclonal antibody to CD274
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr
Immunogen	Full length human CD274
Target	CD274
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (frozen sections): Acetone fixation. Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG2b kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9NZQ7
Entrez	29126
Dilution Range	Immunohistochemistry (frozen sections): Acetone fixation. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of cells stained using anti-h...



Flow cytometry surface staining pattern ...