



Product Datasheet

CD273 antibody (orb402782)


Description	Mouse monoclonal antibody to CD273
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC
Immunogen	human CD273 cDNA
Target	CD273
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG2a kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9BQ51
Entrez	80380
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD11c
positive HLA-
D...