


Product Datasheet

CD49c antibody (orb402783)

Description	Mouse monoclonal antibody to CD49c
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IP, WB
Immunogen	Human SSC-9 cell line (squamous cell carcinoma)
Target	CD49c
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P26006
Entrez	3675
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of
HeLa cells
stained using
a...