



Product Datasheet

Tissue factor antibody (orb402784)


Description	Mouse monoclonal antibody to Tissue factor
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Human brain tissue factor (CD142)
Target	CD142
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 3-12 µg/ml. Western blotting: Recommended dilution: 1-2 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P13726
Entrez	2152
Dilution Range	Flow cytometry: Recommended dilution: 3-12 µg/ml. Western blotting: Recommended dilution: 1-2 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD142
positive
lymph...



Flow
cytometry
surface
staining
pattern ...