



Product Datasheet

TC21 antibody (orb402786)

Description	Mouse monoclonal antibody to TC21
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, WB
Immunogen	Human TC21
Target	R-RAS2
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P62070
Entrez	22800
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Expiration Date	12 months from date of receipt.



Flow cytometry analysis (intracellular staining) of R-Ras2 / TC21 (isoform a) PUR / GAM APC cells.



Immunocytochemistry analysis of colocalization of R-Ras2 / TC21 (isoform a) PUR / GAM APC and GFP in cells.