



Product Datasheet

CD300e antibody (orb402788)

Description	Mouse monoclonal antibody to CD300e
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, IP, WB
Immunogen	CD300e-HA-transfected cells
Target	CD300e
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q496F6
Entrez	342510
Dilution Range	Functional application: Stimulation. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human monocytes (red-fille...



Flow cytometry surface staining pattern ...