



Product Datasheet

CD89 antibody (FITC) (orb402791)

Description

Mouse monoclonal antibody against CD89 conjugated to



Flow
cytometry
analysis
(surface
stainin...

Reactivity

Human, Primate

Conjugation

FITC

Tested

FC

Applications

Immunogen

Ag8.653 myeloma cells

Target

CD89

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.

Isotype

Mouse IgG1 kappa

Clonality

Monoclonal

Purity

Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

[P24071](#)

Entrez

[2204](#)

Dilution Range

Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.

Expiration Date

12 months from date of receipt.