



Product Datasheet

PD1 antibody (orb402797)

Description	Mouse monoclonal antibody to PD1
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr
Immunogen	recombinant human CD279
Target	CD279
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q15116
Entrez	5133
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD279 positive lymph...



Flow cytometry surface staining pattern ...