

Product Datasheet

CD81 antibody (biotin) (orb402798)

Description	Mouse monoclonal antibody against CD81 conjugated to biotin
Reactivity	Feline, Human, Rabbit
Conjugation	Biotin
Tested Applications	FC, ICC, IHC-P, IP, WB
Immunogen	MOLT-4 (human T-ALL cell line)
Target	CD81
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: recommended dilution: 1 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P60033
Entrez	975
Dilution Range	Flow cytometry: recommended dilution: 1 µg/ml.
Expiration Date	12 months from date of receipt.

