


Product Datasheet

PCNA antibody (PE) (orb402809)

Description

Mouse monoclonal antibody against PCNA conjugated to



Separation of
K562 cells
stained using
a...

Reactivity

Drosophila, Gallus, Human, Mouse, Primate, Rat

Conjugation

PE

Tested

FC

Applications

Immunogen

recombinant rat PCNA

Target

PCNA

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Intracellular staining.

Isotype

Mouse IgG2a kappa

Clonality

Monoclonal

Purity

Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

[P12004](#)

Entrez

[5111](#)

Dilution Range

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Intracellular staining.

Expiration Date

12 months from date of receipt.