


Product Datasheet

CD165 antibody (orb402810)

Description	Mouse monoclonal antibody to CD165
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, IP, WB
Immunogen	CD165 purified from human Molt-4 cell line
Target	CD165
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (frozen sections): Application note: acetone fixation. Flow cytometry: Recommended dilution: 1-4 µg/ml. Western blotting: Recommended dilution: 1-2 µg/ml; non-reducing conditions.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Entrez	23449
Dilution Range	Immunohistochemistry (frozen sections): Application note: acetone fixation. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of Jurkat cells stained using...