


Product Datasheet

CD222 antibody (biotin) (orb402816)


Description

Mouse monoclonal antibody against CD222 conjugated

Reactivity	Human, Primate
Conjugation	Biotin
Tested Applications	FC, IP, WB
Immunogen	Recombinant Vaccinia virus encoding CD222.
Target	CD222
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Extracellular and intracellular staining; recommended dilution: 2-6 µg/ml. Western blotting: Recommended dilution: 1 µg/ml; non-reducing conditions.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P11717
Entrez	3482
Dilution Range	Flow cytometry: Extracellular and intracellular staining; recommended dilution: 2-6 µg/ml. Western blotting: Recommended dilution: 1 µg/ml; non-reducing conditions.
Expiration Date	12 months from date of receipt.



Flow cytometry analysis (surface stainin...



Anti-Hu CD222 Biotin (clone MEM-238) wor...