



Product Datasheet

CD222 antibody (orb402817)


Description	Mouse monoclonal antibody to CD222
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	Recombinant Vaccinia virus encoding CD222.
Target	CD222
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Extracellular and intracellular staining. Recommended dilution: 2-6 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P11717
Entrez	3482
Dilution Range	Flow cytometry: Extracellular and intracellular staining. Recommended dilution: 2-6 µg/ml. Western blotting: non-reducing conditions. Recommended dilution: 1-2 µg/ml.
Expiration Date	12 months from date of receipt.



Flow cytometry detection of CD222 in hum...



Separation of human neutrophil granulocy...



Flow cytometry surface staining pattern ...