

Product Datasheet

Human TNFSF11 / RANKL / CD254 Protein (orb402832)

Description

Human TNFSF11 Protein, His Tag is

(HEK293). It contains AA Gly 64 - A...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TNFSF11

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

N-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 22.4 kDa. The protein migrates as 27-35 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

AAC51762.1

Purity

90%

Source

Human TNFSF11, His Tag (active trimer) (orb402832) is expressed from human 293 cells (HEK293). It contains AA Gly 64 - Asp 245 (Accession # AAC51762.1).

MW

22.4 kDa

NCBI

[AAC51762.1](#)

Expiration Date

6 months from date of receipt.