




Product Datasheet

Chat antibody (orb40291)

Description	Rabbit polyclonal antibody to Chat
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	CHAT
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	82536
Uniprot ID	P28329
NCBI	NP_066266.3 , NP_065574.3 , NP_001136406.1 , NP_066264.3 , NP_066265.3 , NP_001136401.1 , NP_001136405.1
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of in NCI-H460 cel...



Western blot analysis of in mouse liver ...