

Product Datasheet

Human LAIR1 / CD305 Protein (orb402999)

Description

Human LAIR-1, Fc Tag is expressed from

Gln 22 - His 163 (Ac...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

LAIR-1

Preservatives

Tris with Glycine, Arginine and NaCl,
pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 42.1 kDa. The protein migrates as 55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

AAH27899.1

Purity

95%

Source

Human LAIR-1, Fc Tag (orb402999) is expressed from human 293 cells (HEK293). It contains AA Gln 22 - His 163 (Accession # Q6GTX8-1).

MW

42.1 kDa

Uniprot ID

[Q6GTX8](#)

NCBI

[AAH27899.1](#)

Expiration Date

6 months from date of receipt.