

Product Datasheet

PRKCB antibody (orb403010)

Description	Goat polyclonal antibody to PRKCB
Species/Host	Goat
Reactivity	Canine, Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PRKCB, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Approx 80kDa band observed in lysates of Human Placenta and Rat Brain (calculated MW of 76.9kDa according to NP_997700.1). Recommended concentration: 0.1-0.3WB: Approx 80kDa band observed in lysates of Human Placenta and Rat Brain (calculated MW of 76.9kDa according to NP_997700.1). Recommended concentration: 0.1-0.3
Protein Sequence	DKEFTRQPVLT.
Clonality	Polyclonal
MW	76.9; 77.0
Entrez	5579
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human placenta ...