


Product Datasheet

C22orf28 antibody (orb403016)

Description	Goat polyclonal antibody to C22orf28
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	HSPC117 (aa201-215), Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Approx 55kDa band observed in lysates of cell lines HeLa, Jurkat and K562 (calculated MW of 55.2kDa according to NP_055121.1). Recommended concentration: 0.1-0.3WB: Approx 55kDa band observed in lysates of cell lines HeLa, Jurkat and K562 (calculated MW of 55.2kDa according to NP_055121.1). Recommended concentration: 0.1-0.3
Protein Sequence	QADPNKVSARAKKR.
Clonality	Polyclonal
MW	55.2
Entrez	51493
Dilution Range	ELISA: 1:64000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using C22orf28 antibody.