


Product Datasheet

ALDH9A1 antibody (orb403020)

Description	Goat polyclonal antibody to ALDH9A1
Species/Host	Goat
Reactivity	Bovine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ALDH9A1, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Approx 50kDa band observed in Human Kidney and Liver lysates (calculated MW of 56.3kDa according to NP_000687.3). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources.WB: Approx 50kDa band observed in Human Kidney and Liver lysates (calculated MW of 56.3kDa according to NP_000687.3). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentrat
Protein Sequence	QKEILDKFTEEVVKQ.
Clonality	Polyclonal
MW	56.3
Entrez	223
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Kidney ly...