

Product Datasheet

Bovine alpha 2 Macroglobulin ELISA Kit (orb403037)

Description	Bovine Alpha-2-macroglobulin(A2M) ELISA kit
Reactivity	Bovine
Range	0.312 µg/mL- 20 µg/mL
Target	A2M
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q7SIH1
Sensitivity	0.078 µg/mL
Expiration Date	Please enquire.