

Product Datasheet

Bovine MSTN ELISA Kit (orb403086)

Description	Bovine Growth/differentiation factor 8(MSTN) ELISA kit
Reactivity	Bovine
Range	6.25 ng/mL-400 ng/mL
Target	MSTN
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O18836
Sensitivity	1.56 ng/mL
Expiration Date	Please enquire.