

Product Datasheet

Bovine HSP90AA1 ELISA Kit (orb403088)

Description	Bovine Heat shock protein HSP 90-alpha(HSP90AA1) ELISA kit
Reactivity	Bovine
Range	6.25 ng/mL-400 ng/mL
Target	HSP90AA1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q76LV2
Sensitivity	1.56 ng/mL
Expiration Date	Please enquire.