

## Product Datasheet

### Bovine PRNP ELISA Kit (orb403115)

<b>Description</b>	Bovine Major prion protein(PRNP) ELISA kit
<b>Reactivity</b>	Bovine
<b>Range</b>	0.01 ng/ml-10 ng/ml
<b>Target</b>	PRNP
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<b>P10279</b>
<b>Sensitivity</b>	0.01 ng/ml
<b>Expiration Date</b>	Please enquire.