

Product Datasheet

Bovine PTX3 ELISA Kit (orb403134)

Description	Bovine Pentraxin-related protein PTX3(PTX3) ELISA kit
Reactivity	Bovine
Range	2.344 ng/mL-150 ng/mL
Target	PTX3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q0VCG9
Sensitivity	0.586 ng/mL
Expiration Date	Please enquire.