

## Product Datasheet

### Bovine VEGFA ELISA Kit (orb403174)

|                          |                                                              |
|--------------------------|--------------------------------------------------------------|
| <b>Description</b>       | Bovine Vascular endothelial growth factor A(VEGFA) ELISA kit |
| <b>Reactivity</b>        | Bovine                                                       |
| <b>Range</b>             | 47 pg/mL-3000 pg/mL                                          |
| <b>Target</b>            | VEGFA                                                        |
| <b>Note</b>              | For research use only                                        |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                 |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                            |
| <b>Uniprot ID</b>        | <b>P15691</b>                                                |
| <b>Sensitivity</b>       | 11.75 pg/mL                                                  |
| <b>Expiration Date</b>   | Please enquire.                                              |