

Product Datasheet

Bovine ORM1 ELISA Kit (orb403191)

Description	Bovine α 1-Acid glycoprotein, α 1-AGP ELISA Kit
Reactivity	Bovine
Range	62.5 ng/mL-4000 ng/mL
Target	ORM1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q3SZR3
Sensitivity	15.6 ng/mL
Expiration Date	Please enquire.