

## Product Datasheet

### Bovine IFN beta ELISA Kit (orb403195)

|                          |                                                             |
|--------------------------|-------------------------------------------------------------|
| <b>Description</b>       | Bovine Interferon $\beta$ , IFN- $\beta$ /IFNB<br>ELISA Kit |
| <b>Reactivity</b>        | Bovine                                                      |
| <b>Range</b>             | Request Information                                         |
| <b>Target</b>            | IFNB1                                                       |
| <b>Note</b>              | For research use only                                       |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                |
| <b>Uniprot ID</b>        | <b>P01578</b>                                               |
| <b>Sensitivity</b>       | Request Information                                         |
| <b>Expiration Date</b>   | Please enquire.                                             |