

Product Datasheet

Bovine MCP1 ELISA Kit (orb403227)

Description	Bovine monocyte chemotactic protein 1/monocyte chemotactic and activating factor,(MCP-1/CCL2/MCAF)EL...
Reactivity	Bovine
Range	78 pg/mL-5000 pg/mL
Target	CCL2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P28291
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.