

Product Datasheet

Bovine SOD1 ELISA Kit (orb403289)

Description	Bovine Cu/Zn-Superoxide Dismutase,Cu/Zn-SOD ELISA Kit
Reactivity	Bovine
Range	15.6 ng/mL-1000 ng/mL
Target	SOD1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma
Uniprot ID	P00442
Sensitivity	3.9 ng/mL
Expiration Date	Please enquire.