

## Product Datasheet

### Bovine SEPP1 ELISA Kit (orb403292)

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Description</b>       | Bovine Selenoprotein P(SEPP1)<br>ELISA kit |
| <b>Reactivity</b>        | Bovine                                     |
| <b>Range</b>             | 10 ng/mL-160 ng/mL                         |
| <b>Target</b>            | SEPP1                                      |
| <b>Note</b>              | For research use only                      |
| <b>Application notes</b> | Detection Wavelength: 450 nm               |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates       |
| <b>Uniprot ID</b>        | <b>P49907</b>                              |
| <b>Sensitivity</b>       | 5 ng/mL                                    |
| <b>Expiration Date</b>   | Please enquire.                            |