

Product Datasheet

Dog HMGB1 ELISA Kit (orb403465)

Description	Dog High mobility group protein B1(HMGB1) ELISA kit
Reactivity	Canine
Range	0.312 ng/mL-20 ng/mL
Target	HMGB1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q6YKA4
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.